

Rotor-Gene® Q accessories

Accessories that can be used with the Rotor-Gene® Q are listed below.

Product	Contents	Cat. no.
36-Well Rotor 	For holding PCR Tubes, 0.2 ml; requires 36-Well Rotor Locking Ring	9018907
36-Well Rotor Locking Ring 	For locking PCR Tubes, 0.2 ml, in the 36-Well Rotor	9018906
72-Well Rotor 	For holding Strip Tubes and Caps, 0.1 ml; requires Locking Ring 72-Well Rotor	9018903
Locking Ring 72-Well Rotor 	For locking Strip Tubes and Caps, 0.1 ml in the 72-Well Rotor	9018904
Rotor-Disc® 72 Rotor 	For holding Rotor-Disc 72 discs in the Rotor-Gene Q; requires Rotor-Disc 72 Locking Ring	9018899



Rotor-Disc 72 Locking Ring	For locking a Rotor-Disc 72 in the Rotor-Disc 72 Rotor	9018900
		
Rotor-Disc 100 Rotor	For holding Rotor-Disc 100 discs in the Rotor-Gene Q; requires Rotor-Disc 100 Locking Ring	9018895
		
Rotor-Disc 100 Locking Ring	For locking a Rotor-Disc 100 in the Rotor-Disc 100 Rotor	9018896
		
Rotor-Disc Pipetting Aid	Aid for marking well during manual reaction setup on a Rotor-Disc Loading Block	9018897
		
Rotor-Disc 72 Starter Kit (115V)	3 Rotor-Disc 72 packs, Rotor-Disc Heat Sealer, Rotor-Disc Heat Sealing Film, Rotor-Disc 72 Rotor and Locking Ring, Rotor-Disc 72 Loading Block, Rotor-Disc Pipetting Aid	9001600
Rotor-Disc 72 Starter Kit (230V)	3 Rotor-Disc 72 packs, Rotor-Disc Heat Sealer, Rotor-Disc Heat Sealing Film, Rotor-Disc 72 Rotor and Locking Ring, Rotor-Disc 72 Loading Block, Rotor-Disc Pipetting Aid	9001601
Rotor-Disc 100 Starter Kit (115V)	2 Rotor-Disc 100 packs, Rotor-Disc Heat Sealer, Rotor-Disc Heat Sealing Film, Rotor-Disc 100 Rotor and Locking Ring, Rotor-Disc 100 Loading Block, Rotor-Disc Pipetting Aid	9001602
Rotor-Disc 100 Starter Kit (230V)	2 Rotor-Disc 100 packs, Rotor-Disc Heat Sealer, Rotor-Disc Heat Sealing Film, Rotor-Disc 100 Rotor and Locking Ring, Rotor-Disc 100 Loading Block, Rotor-Disc Pipetting Aid	9001603
Rotor-Disc Heat Sealer 110V	Heat sealing instrument for use with Rotor-Discs; requires Rotor-Disc 72 or 100 Loading Block	9018898
		

Rotor-Disc Heat Sealer 230V



Heat sealing instrument for use with Rotor-Discs; requires Rotor-Disc 72 or 100 Loading Block

9019725

Loading Block 72 x 0.1 ml Tubes



Aluminum block for manual reaction setup with a single-channel pipet in 72 x 0.1 ml tubes


9018901

Loading Block 72 x 0.1 ml Multi-channel



Aluminum block for reaction set-up with multichannel pipets in 72 x 0.1 ml tubes 9018902

Loading Block 96 x 0.2 ml Tubes



Aluminum block for manual reaction set-up in a standard 8 x 12 array using 96 x 0.2 ml tubes

9018905


Rotor Holder



Metal free-standing holder for assembling tubes and Rotor-Discs into rotors

9018908


Rotor-Disc 72 Loading Block



Aluminum block for manual and automated reaction setup in Rotor-Disc 72 discs

9018910

Rotor-Disc 100 Loading Block



Aluminum block for manual and automated reaction setup
in Rotor-Disc 100 discs

9018909

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, (QIAGEN Group); Rotor-Gene®, Rotor-Disc™ (Corbett Research Pty Ltd).

06/2009 © 2009 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

China ■ 021-3865-3865

Denmark ■ 80-885945

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Austria ■ 0800/281010

France ■ 01-60-920-930

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 1544 7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-639

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 65-67775366

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample & Assay Technologies